

# QCPP Compliant Dispensary Refrigerators



Quality Care  
Pharmacy Program

QCPP recognises purpose-built dispensary refrigerators on this list as meeting the requirement of 'Compliant Dispensary Refrigerator' in Element 5 Action 4.

QCPP acknowledges that refrigerators that do not appear on this list may still comply with the QCPP requirements. If a pharmacy receives a non-conformance at QCPP Assessment because the dispensary refrigerator(s) is not on this list, we encourage submission of a Remedial Review Request at [www.qcpp.com/resources/tools/qcpp-forms](http://www.qcpp.com/resources/tools/qcpp-forms).

Company		Compliant models	
<b>All-Lab Scientific</b>	ALV-358		
<b>Arcus Australia Pty Ltd</b> 08 9381 3888 <a href="http://www.arcus-au.com.au">www.arcus-au.com.au</a>	SV PV	AUD1-SV AUD1-V CMU	
<b>Arrowsmith and Grant Refrigeration</b> 03 9706 6655 <a href="http://www.arrowsmithgrant.com.au">www.arrowsmithgrant.com.au</a>	AG 220 V		
<b>Avem Quirks Medisafe Refrigerator</b> 1800 658 120 <a href="http://www.avemquirks.com.au">www.avemquirks.com.au</a>	<b>Medisafe:</b> 125 311 371 (glass or solid door) 410 471 BUG FM1/1 HB80 R88 (small) S5 TME1000-C FKG 371	<b>Vestfrost:</b> BC85 SLC125	
<b>BD Rowa</b> 1800 833 518 <a href="http://www.carefusion.com.au">www.carefusion.com.au</a>	Rowa Vmax system		
<b>Bromic MediFridge®</b> 1300 276 642 <a href="http://www.bromicrefrigeration.com.au">www.bromicrefrigeration.com.au</a>	MED0140GD MED0374GD MED0220GD	VSS 374 T MediFridge 145L	
<b>B&amp;R Refrigeration Solutions</b> 0430 280 865 <a href="http://www.brrefrigerationsolutions.com">www.brrefrigerationsolutions.com</a>	BRSC 125 BRSC 300 BRSC 400	HYC-L130	

<b>Crown Scientific</b>	<b>Dometic:</b> MP 150 MP 155 MP 285 MP 405		
<b>CSK Sales (Pinnacle)</b> 07 3376 5277 www.cskgroup.com.au	PPR-S-66-1-GD PPR-S-145-1-GD PPR-S-300-1-GD PPR-S-430-1-GD PPR-M-1081-2-GD	CSK PPR 145 CSK PPR 430	
<b>Elk Scientific</b> 0423 520 302 www.elkscientific.com.au	<b>Logicold Pharmacy:</b> ELKLCP30 ELKLCP60 ELKLCP125 ELKLCP300 ELKLCP400		
<b>Euro Chill</b> 1300 459 140 www.eurochill.com.au	Vacc Safe 80 PLUS Vacc Safe 311 PLUS Vacc Safe 371 PLUS Vacc Safe 120 PREMIUM Vacc Safe 250 PREMIUM Vacc Safe 400 PREMIUM Vacc Safe 600 PREMIUM Vacc Safe 1000 PREMIUM Vacc Safe Premium 120 EC	VS80 EC-80 Ltr VS311 EC-310 Ltr VS370 EC-370 Ltr VS250 EC VS400 EC Vacc Safe VS80 EC Vacc Safe VS311 EC Vacc Safe VS371 EC Vacc S Select 311 EC VS Select 370 EC	<b>Mediguard:</b> Mediguard 60 Mediguard 401 Mediguard 601
<b>Exquisite Marketing Australia</b> 03 9561 1595	MV30 MV60 MV150 MV300 MV400		
<b>ICS Pacific</b> 03 9706 5115 www.icspacific.com.au	Pharma 1000GD Pharma 3000GD Pharma 3000SD Pharma 3	G105L Medical G125L Medical G130L Medical G135L Medical G315L Medical G375L Medical G400L Medical G410L Medical	F135 F375P
<b>LabGear Australia</b> 07 3368 1112 www.labgearaustralia.com.au	<b>Pharm Chill:</b> BC-29 BC-140 BC-300 BC-400		

<b>Laboratory Equipment Pty Ltd</b> 02 9560 2811 www.labec.com.au	LABEC PVR-130GD LABEC PVR-440GD	LABEC HR200GP LABEC PVR-374 LABEC GMO374	
<b>Liebherr Mediline (Andi-Co Australia Pty Ltd)</b> 1800 685 899 www.liebherrprofessional.com.au	LKUv1610 LKUv1612 LKUv1613 LKv3910 LKv3912 LKv3913	FKVSL3612	
<b>MediLine Equipment</b> 0400 185 682	125ML 300ML 400ML	<b>Labcold:</b> EC110 EC300 EC400	
<b>Nuline Refrigeration</b> 03 9562 3355 info@nulinerefrigeration.com.au	HR200 HR400 HR600 HR400P HR700P HR400G HR400GP	BEV32P BEV40P BEV60P Mondial Elite BEVPV40 Elite 310	GN1420TNP MF1210TNG MFT1410TNG ML400 NLMB/2 NLMB/2 P NLM700-1P NLM1000/2 P NM1000/2 NMM500/1 NMM1000/2 TTG-14P TTGPR14L
<b>Orford Refrigeration Pty Ltd</b> 07 4639 9000 www.orford.com.au	V V D P07A (small) V V D P09A (med) V V D P14A (large)	FM20S2-R-A FM36S2-R-A FM45S2-R-A	
<b>Quantum Scientific</b>	<b>Sanyo Medi-Cool:</b> MPR-161D(H) MPR-311D(H) MPR-514D (H)	<b>Sanyo Labcool:</b> MPR-721	
<b>Rapini Pty Ltd</b> 03 9555 0848 www.rapini.com.au	R155-G R300-G	WMM1000/2P	

<p><b>Rollex Vaccine Refrigerators</b>          Rollex Group (Australia) Pty Ltd          1300 880 441          www.rollexmedical.com.au</p>	<p><b>LEC:</b>          350          503          KS290          PE102          PE109          PE202          PE207          PE350          PE357          PE501          PE502          PE507          PE900          PE 901          PE902          PE907          PE1202          PE1602          PE1607          PG102          PG202          PG207          PG357          PG 501          PG502          PG507          PG901          PG902          PG907          PG1202          PG1602          PG1607          PGR151          PGR273          PGR334          PGR353          PSR151          PSR273          PSR353          PSR334          PT555</p>	<p><b>Pol Eko - MATOS1*:</b>          MATOS-eco 68 R/G          MATOS-eco 150 R          MATOS-eco 150 R/G          MATOS-eco 200 R          MATOS-eco 200 R/G          MATOS-eco 250 R          MATOS-eco 300 R          MATOS-eco 300 R/G          MATOS-eco 493 R          MATOS-eco 493 R/G          MATOS-eco 625 R          MATOS-eco 625 R/G          MATOS-eco 1365 R/DGT          MATOS-eco 1365 R/GDT          MATOS-eco 1460 R/G          MATOS-eco PLUS 250          glass/solid door          MATOS-eco PLUS Cloud          400 glass/solid door R/G          MATOS-eco PLUS Cloud          1365 R/DT          MATOS-eco PLUS Cloud          1460 R/DT &amp; R/GDT          CHL 1+ A          CHL 2+          CHL 2+ A          CHL 3+ A          CHL 4+          CHL 5+          CHL 5+ A          CHL 500          CHL 500 A          CHL 700          CHL 700 A          CHL 1200          CHL 1450 TROPIC</p> <p><b>Aria – MATOS</b>          MATOS Aria-Eco Cloud          150 glass/solid door          MATOS Aria-Eco 220          glass door/solid door          MATOS Aria-Eco 300          glass door/solid door          MATSO Aria-Eco 300          Cloud solid door          MATOS Aria-Eco 374          glass door/solid door          MATOS Aria Eco 374          Cloud solid</p>	<p><b>RK - Koldpia:</b>          RK600G          RK600S          RK1100G          RK1100S          SGRR-23FD          SGRR-23GD          SGRR-23SD          SGRR-49FD          SGRR-49GD          SGRR-49SD</p> <p><b>Angelantoni:</b>          FRL 360          FRL 360 Glass          FRL 500          FRL 500 Glass          700 TN          700 TN-GL          1500 TN          1500 TN-GL</p> <p><b>Twinbird Portable:</b>          SC-BV25</p>
--	---	---	--

<sup>1</sup> \*MATOS-cloud models are a variation of MATOS-eco, which include online monitoring capability and are also compliant.

<b>Skope Industries</b>	TME475XL-CK TME100-H TME650-H TME1000-C TME1500-C	B130-PH B170E-PH B310E-PH RCTME1000	LM620/665/190 SK650-2-CK SK1300-2-C SK1500XLP-C
<b>Thermoline Scientific Refrigerators</b> 02 9604 3911 www.thermoline.com.au Note: 1-GD = 1 glass door 2-GD = 2 glass doors 3-GD = 3 glass doors	TEPR-60-1-GD TEPR-80-1-GD TEPR-140-1-GD TEPR-145-1-GD TEPR-200-1-GD TEPR-230-1-GD TEPR-250-1-GD TEPR-400-1-GD TEPR-440-1-GD TEPR-500-1-GD TEPR-740-2-GD TEPR-142S-R	TLMPR140-1-GD TLMPR220-1-GD TLMPR300-1-GD TLMPR350-1-GD TLMPR400-1-GD TLMPR800-2-GD TLMPR900-2-GD TLMPR1215-3-GD TLMPR1275-3-GD	TPR-250-1-GD TPR-520-1-GD TPR-530-1-GD TPR-750-2-GD TPR-800-2-GD TPR-900-2-GD TPR-950-2-GD TPR1086-2-GD
<b>Vaccine Fridges Australia (Frigidum Pty Ltd)</b> 1300 272 464 www.vaccinefridges.com.au	Huxford VXH20		
<b>VWR International Pty Ltd</b> 1300 727 696 csaus@au.vwr.com	<b>Panasonic:</b> MPR-S163-PE	<b>Dometic:</b> MP 155 MP 285 MP 405	<b>Sanyo Medi-Cool:</b> MPR-161D (H) MPR-311D (H) MPR-514D (H)
<b>Yili Ramli &amp; Co. Pty Ltd</b> 02 9679 9488 www.vaccinerefrigeration.com.au	VaccineCab YC-280		
<b>Other</b>	<b>Linder and May</b> LMPR-400-1-GD  <b>Tropical Thermaster</b> LG-570GTH TG-570GTH	<b>Baker Refrigeration</b> Custom built  <b>Norcool</b> Super 5	<b>Rowsen</b> BYY-300